

Cobalt(II) protein NMR

*This is a project currently being explored by myself between Chemistry and Agriculture
The protein is a cobalt substituted peptide deformylase, an enzyme required for the
growth of bacteria and weed. We are trying to find an effective inhibitor to this enzyme as
antibiotics and herbicides.*

Step 1. Record a normal 1D spectrum of Co_DEF on a 400 MHz machine.

Setup the experiment using s2pul sequence;

nt=1, tpwr=63

ga

aph

Place the left cursor on the water peak, **nl**

movetof

ga


array

pw *array pw90 and set pw=pw90.*

double check *tof* by typing **movetof**

ga

you should see the following spectrum.



The parameters used for collecting this spectrum are listed as follows:

ACQUISITION		SAMPLE		PROCESSING			FLAGS	
sfrq	399.730	date	Jun 1 2004	lb	not used	il		n
tn	H1	solvent	D2O	sb	not used	in		n
at	3.744	file	exp	gf	not used	dp		y
np	44932	DECOUPLING		awc	not used	hs		nn
sw	6000.6	dn	H1	lsfid	not used	temp	SPECIAL	not used
fb	3000	dof	0	phfid	not used			
bs	16	dm	nnn	wtfile				
ss	0	dmm	c	proc		ft		
tpwr	60	dmf	11148	fn	not used			
pw	16.2	dpwr	42	math		f		
p1	0							
d1	1.000			werr				
d2	0			wexp				
tof	-153.9			wbs				
nt	1			wnt				
ct	0							

Please write down *tof* which is *-153.9* and *pw* which is *16.2*.

Step2. Use WEFT sequence to run a protein NMR.

jexp2 ↵

we have a sequence called weftafhop, you are supposed to have one by contact you manager.

Load the original fid of a sample someone previously used, or the parameters.

Insert your sample(90%H2O and 10%D2O).

dg ↵

tpwr=63 nt=2048 ↵


pw90=11.5 ↵ (*tpwr* going up 60→63; so *pw90* going down 16.2/1.414)

(If you want to know why, go to the *pw90* calibration chapter)

pw=pw90 ↵

supfreq=tof ↵

dps ↵




dg ↵

ACQUISITION	SAMPLE	PROCESSING	FLAGS
sfrq 399.730	seqfil weftafhop	lb not used	il n
tn H1	date Jun 2 2004	sb not used	in n
at 0.025	sample WT Fe50Dred	gf 0.003	dp y
np 6250	solvent D20	gfs not used	pregrad n
sw 125000.0	file exp	awc not used	prex n
fb 69000	DECOUPLING	lsfid -2	shape n
bs 16	dn H1	phfid not used	weftpls y
ss 256	dof 26321.5	proc lp	gradflg y
tpwr 63	dm nnn	fn not used	satflg y
pw 11.5	dmm c		satpwr 8
p1 23.6	dmf 16949	shape n	gt0 0.001000
d1 0.010	dpwr 8	shapwr 0	gzlv10 5000
d2 0	SPECIAL	supfreq -153.9	pp 8.05
tof -153.9	temp not used	pp 8.05	gt1 0.000200
nt 65536	gain 30	pw90 11	weftdly 0.035000
ct 0			gradly 0.000100

Things you need to take care of are:
gain not exceed 30, *satpwr* smaller than 12.


ga ↵



After the acquisition is done, you are supposed to have the above spectrum.

The inhibitor would react with the metal and stop the functioning of the protein. Here actinolin is used to chelate the cobalt atom. The same parameters are used for the NMR study.

A spectrum is shown on the next page.



To do this job more informatively, you may do the following to optimize the parameters.

1. **array supfreq (in the range ± 15 is enough, because it is well set in the normal 1D)**
2. **array (rof2,alfa)**
3. **array (d1,weftdly)**

Keep in mind that **weftdly = at + d1**

You should manually phase the spectrum. Keep practicing and you will be able to phase it. Phasing error may occur.

To check the error, type **lp?**; if lp is far away from zero, then you need to correct it by typing **crof2**
ga

You will get a better spectrum.

- **rof1**: receiver gating time preceding pulse 0-8190;
- **rof2**: receiver gating time following pulse 0-8190;
- **alfa**: After the final event in the pulse sequence, including any receiver gate times occurring following the pulse, acquisition occurs after a delay. This delay includes a fixed part, **alfa**, and a variable part, **1/(beta*fb)**.

How to array rof2 and alfa

```
rof2? ↵ ( rof2=1)
alfa? ↵ (alfa=7)
array='(rof2,alfa)'
rof2=0,1,2,3,4,5,6,7,8
alfa=8,7,6,5,4,3,2,1,0
```

(remember: rof2 + alfa = 8 all the time)

How to array d1 and weftdly

```
at? ↵ or dg↵ to see what it is. at=0.025
array='(d1,weft)' ↵
d1=0.01,0.03,0.06,0.1↵
weft=0.035,0.055,0.085,0.125↵
```

(see that at + d1 =weftdly)

Hint:

1. Do not miss the brackets, “(“ and “)” when you type array='(xxx,yyy)' in this case, because if they are missing, you are arraying in 9×9 dimension which will run 81 jobs. You only need 9 jobs. Please type **da** to see the details in the text. **ds(1)** to display the first one, **ds(4)** to display the 4th spectrum. Compare all the spectra and select the best one. Use the best parameter to run protein NMR.
2. Use pure D2O as solvent to assign the peaks to the histidine moiety. To do this, you need to change the solvent and record another spectrum.